

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-DHRS-403**


# Anti DHRS

|                           |                                                                                                                                                                                                      |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product type</b>       | Primary antibodies                                                                                                                                                                                   |
| <b>Immunogen</b>          | DHRS11 (ENSP00000354182, 1 a.a. ~ 181 a.a) full-length human protein.                                                                                                                                |
| <b>Form</b>               | Liquid<br>Protein A purified                                                                                                                                                                         |
| <b>Volume</b>             | 100ul                                                                                                                                                                                                |
| <b>Specificity</b>        | Human                                                                                                                                                                                                |
| <b>Host</b>               | Rabbit                                                                                                                                                                                               |
| <b>Immunogen Sequence</b> | MFSAIRSQHSGVDICINNAGLARPD TLLSGSTSGWKDMFNVNVLALSICTREAYQSMKER<br>NVDDGHIININSMSGHRVLP LSVTHFY SATKYAVTALTEGLRQELREAQTHIRATCISPG<br>VVETQFAFKLHDKDPEKAAATYEQMKCLKPEDVAEAVIYVLSTPAHIQIGDIQMRPTEQV<br>T |

ARPG836, FLJ39232, MGC4172, SDR24C1

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~20.00 kDa using DHRS11 transfected 293T lysate.

**Western Blotting (Transfected lysate)**




Western Blot analysis of DHRS11 expression in transfected 293T cell line by DHRS11 rabbit polyclonal antibody.

Lane 1: DHRS11 transfected lysate(20.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DHRS11 expression in mouse intestine (Predicted Mw: 28.3080 kDa).

**Western Blotting (Tissue lysate)**



DHRS11 rabbit polyclonal antibody. Western Blot analysis of DHRS11 expression in mouse intestine.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619