

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-DLAT-028**

# Anti DLAT

|                     |                                                                  |
|---------------------|------------------------------------------------------------------|
| <b>Product type</b> | Primary antibodies                                               |
| <b>Immunogen</b>    | DLAT (NP_001922.2, 1 a.a. ~ 647 a.a) full-length human protein.. |
| <b>Form</b>         | Liquid                                                           |
|                     | Protein A purified                                               |
| <b>Volume</b>       | 100ul                                                            |
| <b>Specificity</b>  | Human                                                            |
| <b>Host</b>         | Rabbit                                                           |


**Immunogen Sequence**

MWRVCARRAQNVAPWAGLEARWTALQEVPGTPRVTSRSGPAPARRNSVTTGYGGVRLCG  
WTPSSGATPRNRLLLQLLGSPGRRYSLPPHQKVLPLSLPTMQAGTIARWEKKEGDKIN  
EGDLIAEVETDKATVGFESLEECYMAKILVAEGTRDVPIGAIICITVGKPEDIEAFKNYT  
LDSSAAPTQQAAPAPTPAATASPPTPSAQAPGSSYPHMQVLLPALSPMTMTMGTVQRWEK  
KVGEKLSEGDLLAEIETDKATIGFEVQEEGYLAKILVPEGTRDVPLGTPLCIIVEKEADI  
SAFADYRPTEVTDLKPQVPPPTPPPVAAPVPTPQPLAPTPSAPCPATPAGPKGRVFSPL  
AKKLAVEKGIDLTQVKGTGPDGRITKKDIDSFVPSKVAPAPAAVVPPTGPGMAPVPTGVF  
TDIPISNIRRVIAQRLMQSKQTIPHYLLSIDVNMGEVLLVRKELNKILEGRSKISVNDFI  
IKASALACKLVPEANSSWMDTVIRQNHVVDVSVAVSTPAGLITPIVFNAHIKGVETIAND  
VVSLATKAREGKLQPHEFQGGTFTISNLGMFGIKNFSAINPPQACILAIGASEDKLVPA  
DNEKGFDVASMMSVTLSCDHRVVDGAVGAQWLAEFRKYLEKPITMLL

**Gene Alias** DLTA, PDC-E2, PDCE2

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~69.00 kDa using DLAT transfected 293T lysate.

**Western Blotting (Transfected lysate)**




Western Blot analysis of DLAT expression in transfected 293T cell line by DLAT rabbit polyclonal antibody.

Lane 1: DLAT transfected lysate(69.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DLAT expression in A-431 (Predicted Mw: 65.7810 kDa).

**Western Blotting (Cell lysate)**



DLAT rabbit polyclonal antibody. Western Blot analysis of DLAT expression in A-431.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619