

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-ESF1-366

Anti ESF1

Product type	Primary antibodies
Immunogen	ESF1 (AAI46543.1, 1 a.a. ~ 851 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSSKQEIMSDQRFRRVAKDPRFWEMPEKDRKVKIDKRFRAMFHDKKFKLNYAVDKRGRPI SHSTTEDLKRFDLSDSDSNLSGEDSKALSQKKIKKKKTQTKKEIDSKNLVEKKKETKKA NHKGSENKTDLDNSIGIKKMKTSCKFKIDSNI SPKKDSKEFTQKNKKEKNIVQHTTDS LEEKQRTLDSGTSEIVKSPRIECSKTRREM QS VVQLIMTRDSDGYENSTDGEMCDKDALE EDSESVSEIGSDEESEN EITSVGRASGDDDGSEDEDEDEDEDEDEDEDEDEDEDEDEDEDE DLARGKGN IETSS EDEDDTADLFPEESGFEHAWRELDKDAPRADEITRRLAVCNMDWDRL KAKDLLALFNSFKPKGGVIFSVKIYPSEFGKERMKEEQVQGPVELLSIPEDAPEKDWT SR EKLRDYQFKRLKYYYAVVDCDSPETASKIYEDCDGLEFESSCSFIDLRFPDDITFDDEP KDVASEVNLTAYKPKYFTSAAMGTSTVEITWDETDHERITMLNRKFKKEELLDMDFQAYL ASSSEDEEEIEEELQGDDGVNVEEDGKTKKSQKDDEEQIAKYRQLLQVIQEKEKKGKEND MEMEIKWV PGLKESAEEMVKNKLEGKDKLTPWEQFLEKKKEKKRLKRKQKALAE EASEEE LPSDVDLNDPYFAEEVKQIGINKKSVKSAKDGTSP EEEIEIERQKAEMALLMDEDED SK KHFNYNKIVEHQNL SKKKKKQLMKKKELIEDDFEVNVNDARFQAMYTSHLFNLDPSDPNF KKT KAMEKILEEKARQ RERKEQELTQA IKKK ESEIEKESQRKSIDPAL SMLIKSIKTKTE QFQARKKQKVK


ABTAP, C20orf6, FLJ20368, HDCMC28P, bA526K24.1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~93.61 kDa using ESF1 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of ESF1 expression in transfected 293T cell line by ESF1 rabbit polyclonal antibody.

Lane 1: ESF1 transfected lysate(93.61 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619