

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-FABP7-358

Anti FABP7

Product type	Primary antibodies
Immunogen	FABP7 (AAH12299.1, 1 a.a. ~ 132 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MVEAFCATWKLNSQNFDEYMKALGVGFATRQVGNVTKPTVIISQEGDKVVIRTLSTFKN TEISFQLGEEFDETTADDRNCKSVVSLDGDKLVIHQKWDGKETNFVREIKDGKMVMTLTF GDVAVVRHYEKA


B-FABP, BLBP, DKFZp547J2313, FABPB, MRG

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~14.63 kDa using FABP7 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of FABP7 expression in transfected 293T cell line by FABP7 rabbit polyclonal antibody.


Lane 1: FABP7 transfected lysate(14.63 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of FABP7 transfected lysate..

Immunoprecipitation



Immunoprecipitation of FABP7 transfected lysate using FABP7 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with FABP7 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619