



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-FKBP11-206


Anti FKBP11

Product type	Primary antibodies
Immunogen	FKBP11 (NP_057678.1, 1 a.a. ~ 201 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MTLRPSLLPLHLLLLLLLLSAAVCRAEAGLETESPVRTLQVETLVEPPEPCAEPAAFGDTL HIHYTGSGLVDGRIIDTSLTRDPLVIELGQKQVIPGLEQSLDMCVGEKRRRAIIPSHLAYG KRGFPPSVPADAVVQYDVELIALIRANYWLKLVKGILPLVGMAMVPALLGLIGYHLYRKA NRPKVSKKKLKEEKRNSKKK

Gene Alias FKBP19, MGC54182

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~22.20 kDa using FKBP11 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of FKBP11 expression in transfected 293T cell line by FKBP11 rabbit polyclonal antibody.


Lane 1: FKBP11 transfected lysate(22.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of FKBP11 transfected lysate.

Immunoprecipitation



Immunoprecipitation of FKBP11 transfected lysate using FKBP11 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with FKBP11 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619