

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-GALNS-246

Anti GALNS

Product type	Primary antibodies
Immunogen	GALNS (NP_000503.1, 1 a.a. ~ 522 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAAVVAATRWWQLLLVLSAAGMGASGAPQPPNILLLLMDDMGWGLGVYGEPSRETPNLD RMAAEGLLFPNFYSANPLCSPSRAALLTGRLPIRNGFYTTNAHARNAYTPQEIVGGIPDS EQLLPELLKKAGYVSKIVGKWHLGHRPQFHPLKHGFDEWFGSPNCHFGPYDNKARNIPV YRDWEMVGRYEEFPINLKTGEANLTQIYLQEALDFIKRQARHHPFFLYWAVDATHAPVY ASKPFLGTSQRGRYGDVREIDDSIGKILELLQDLHVADNTFVFFTSDNAGAALISAPEQG GSNGPFLCGKQTTFEGGMREPALAWWPGHVTAGQVSHQLGSIMDLFTTSLALAGLTPPSD RAIDGLNLLPTLLQGRMLMDRPIFYRGTLMATLGQHKAHFWTWTNSWENFRQGIDFCP GQNVSGVTTNLEDHTKLPLIFHLGRDPGERFPLSFASA EYQEALSRITSVVQQHQEALV PAQPQLNVCNWAVMNWAPPGCEKLGKCLTPPESIPKKCLWSH
Gene Alias	FLJ17434, FLJ42844, FLJ98217, GALNAC6S, GAS, MPS4A
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~58.00 kDa using GALNS transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of GALNS expression in transfected 293T cell line by GALNS rabbit polyclonal antibody.

Lane 1: GALNS transfected lysate(58.00 KDa).

Lane 2: Non-transfected lysate.

IP Immunoprecipitation



Immunoprecipitation of GALNS transfected lysate using GALNS rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GALNS rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619