

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GALNT12-374

Anti GALNT12

Product type	Primary antibodies
Immunogen	GALNT12 (AAH13945.1, 1 a.a. ~ 272 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAGGLFAVSKKYFEYLGSYDTGMEVWGGENLEFSFRIWQCGGVLETHPCSHVGHVFPKQA PYSRNKALANSVRAAEVWMDEFKELYHRNPRARLEPFGDVTERKQLRDKLQCKDFKWFL ETVYPELHVPEDRPGFFGMLQNKGLTDYCFDYNPPDENQIVGHQVILYLCHGMGQNQFFE YTSQKEIRYNTHQPEGCIAVEAGMDTLIMHLCEETAPENQKFILQEDGSLFHEQSKKCVQ AARKESDSFVPLLRDCTNSDHQKWFFKERML


FLJ21212, GalNAc-T12

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~31.70 kDa using GALNT12 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of GALNT12 expression in transfected 293T cell line by GALNT12 rabbit polyclonal antibody.

Lane 1: GALNT12 transfected lysate(31.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of GALNT12 expression in Jurkat (Predicted Mw: 66.9380 kDa).

Western Blotting (Cell lysate)



GALNT12 rabbit polyclonal antibody. Western Blot analysis of GALNT12 expression in Jurkat.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619