

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-GCDH-232

Anti GCDH

Product type	Primary antibodies
Immunogen	GCDH (NP_000150.1, 1 a.a. ~ 438 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MALRGVSVRLLSRGPGLHVLRTWVSSAAQTEKGGRTQSQLAKSSRPEFDWQDPLVLEEQL TTDEILIRDTFRTYCQERLMPRILLANRNEVFHREIISEMGELGVLGPTIKGYGCAGVSS VAYGLLARELERVDSGYRSAMSVQSSLVMHPIYAYGSEEQRQKYLPQLAKGELLGCFGLT EPNSGSDPSSMETRAHYNSSNKSYSYTLNGTKTWITNSPMADLFVWWARCEDGCIRGFLLEK GMRGLSAPRIQKGKFSLRASATGMIIMDGVEVPEENVLPGASSLGGPFGCLNNARYGIAWG VLGASEFCLHTARQYALDRMQFGVPLARNQLIQKKLADMLTEITLGLHACLQLGRLKDQD KAAPEMVSLKRNCGKALDIARQARDMLGGNGISDEYHVIRHAMNLEAVNTYEGTHDIH ALILGRAITGIQAFTASK
Gene Alias	ACAD5, GCD
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~48.10 kDa using GCDH transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of GCDH expression in transfected 293T cell line by GCDH rabbit polyclonal antibody.

Lane 1: GCDH transfected lysate(48.10 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of GCDH expression in human kidney (Predicted Mw: 48.1270 kDa).

Western Blotting (Tissue lysate)



GCDH rabbit polyclonal antibody. Western Blot analysis of GCDH expression in human kidney.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619