

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-GCNT1-222

Anti GCNT1

Product type	Primary antibodies
Immunogen	GCNT1 (NP_001481.2, 1 a.a. ~ 428 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MLRTLRRRLFSYPTKYYFMVLVLSLITFSLRIHQKPEFVSVRHLELAGENPSSDINCT KVLQGDVNEIQVKLEILTVKFKRPRWTPDDYINMTSDCSSFIKRRKYIVEPLSKEEAE FPIAYSIVVHHKIEMLDRLLRAIYMPQNFYCIHVDTKSEDSYLAAVMGIASCFSNVFVAS RLESVVYASWSRVQADLNCMKDLYAMSANWKYLINLCGMDFPIKTNLEIVRKLKLLMGEN NLETERMPSHKEERWKKRYEVVNGKLTNTGTVKMLPPLETPLFSGSAYFVVSREYVGYYL QNEKIQKLMEWAQDTYSPDEYLWATIQRIPEVPGSLPASHKYDLSDMQAVARFVKWQYFE GDVSKGAPYPPCDGVHVRVSVCFGAGDLNWMLRKHHLFANKFDVDVDLFAIQCLDEHLRH KALETCLKH
Gene Alias	C2GNT, C2GNT-L, C2GNT1, G6NT, MGC126335, MGC126336, NACGT2, NAGCT2
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.80 kDa using GCNT1 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of GCNT1 expression in transfected 293T cell line by GCNT1 rabbit polyclonal antibody.

Lane 1: GCNT1 transfected lysate(49.80 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of GCNT1 expression in human pancreas (Predicted Mw: 49.7850 kDa).

Western Blotting (Tissue lysate)



GCNT1 rabbit polyclonal antibody. Western Blot analysis of GCNT1 expression in human pancreas.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619