



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GOSR1-015

Anti GOSR1

Product type Primary antibodies

Immunogen GOSR1 (NP_004862.1, 1 a.a. ~ 250 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MAAGTSSYWEDLRKQARQLENELDLKLVSFSLCTSYSHSSTRDGRRDRYSSDTTPLLNG
SSQDRMFETMAIEIEQLLARLTGVNDKMAEYTNAGVPSLNAALMHTLQRHRDILQDYTH
EFHKTKANFMAIRERENLMGSVRKDIESYKSGGVNNRRTELFLKEHDHLRNSDRLEET
ISIAMATKENMTSQRGMLKSIHSMNTLANRFPVNSLIQRINLRKRRDSLILGGVIGIC
TILLLLYAFH

Gene Alias GOLIM2, GOS-28, GOS28, GOS28/P28, GS28, P28

Application notes IP
Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of GOSR1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of GOSR1 transfected lysate using GOSR1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GOSR1 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619