

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GPAM-018

Anti GPAM

Product type	Primary antibodies
Immunogen	GPAM (AAH30783.1, 1 a.a. ~ 828 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

**Immunogen
Sequence**

MLPPQKKPWESMAKGLVLGALFTSFLLLMYSYAVPPLHAGLASTTPEAAASCSPPALEPEMD
EYALT LGTIDVSYLPHSSEYSVGRCKHTSEEWGECGFRPTIFRSATLKWKESLMSRKR
PFVGRCCYSCTPQSWDKFFNTSIPSLGLRNVIIYINETHTRHRGWLARRLSYVLFIQERDV
HKGMFATNVTENVLNSSRVQEAIAEVAEELNPDGSAQQQSKAVNKVKKKAKRILQEMVAT
VSPAMIRLTGWVLLKLFNSFFWNIQIHKGQLEMVKAATETNLP LLFLPVHRSHIDYLLLT
FILFCHNIKAPYIASGNNLNIPFSTLIHKLGGFFIRRRLDETPDGRKDVLYRALLHGH
VELLRQQQFLEIFLEGTRSRSGKTSCARAGLLSVVVDLSTNVIPDILIPVGISYDRII
EGHYNGEQLGKPKKNESLWSVARGVIRMLRKNYGCVRVDFAPFSLKEYLESQSQKPVSA
LLSLEQALLPAILPSRPSDADEGRDTSINESRNATDESLRRRLIANLAEHILFTASKSC
AIMSTHIVACLLLYRHRQGIDLSTLVEDFFVMKEEVLARDFDLGFSGNSEDEVVMHAIQLL
GNCVTITHTSRNDEFFITPSTTVPSVFELNFYSNGVLHVFIMEAIIACSLYAVLNKRGLG
GPTSTPPNLISQEQVLRKAASLCYLLSNEGTLSPCQTFYQVCHETVGKFIQYGILTVAE
HDDQEDISPLAEQQWDKKLPEPLSWRSDEEDEDSDFGEEQRDCYLKVSQSKEHQQFITF
LQRLGPLEAYSSAAIFVHNFSGPVPEPEYLQKLHKYLITRTERNVAVYAESATYCLVK
NAVKMFKDIGVFKETKQKRVSVLELSSTFLPQCNRQKLLEYILSFVVL


Gene Alias GPAT1, KIAA1560, MGC26846, RP11-426E5.2

**Application
notes**

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~93.90 kDa using GPAM transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of GPAM expression in transfected 293T cell line by GPAM rabbit polyclonal antibody.

Lane 1: GPAM transfected lysate(93.90 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619