



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GPN1-542


Anti GPN1

Product type	Primary antibodies
Immunogen	GPN1 (NP_009197.1, 1 a.a. ~ 374 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAASAAAELQASGGPRHPVCLLVLMAGSGKTTFVQRLTGHLHAQGTPPYVINLDPVH EVPFPANIDIRDTVKEYEVMKQYGLGPNGGIVTSLNLFATRFDQVMKFIEKAQNMSKYVL IDTPGQIEVFTWSASGTIITEALASSFPTVVIYVMDTSRSTNPVTFMSNMLYACSIKYT KLPIFVVMNKTDIIDHSFAVEWMQDFEAFQDALNQETTYVSNLTRMSLVLDEFYSSLRV VGVSAVLGTGLDELQVTSAAEEYEREYRPEYERLKKSLANAESQQQREQRLRLRKDMG SVALDAGTAKDSLSPVLHPSDLILTRGTLDEEDEEADSDTDDIDHRVTEESHEEPAFQNF MQESMAQYWKRNNK

Gene Alias ATPBD1A, MBDIN, NTPBP, XAB1

Application notes **IP**
Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of GPN1 transfected lysate

Immunoprecipitation



Immunoprecipitation of GPN1 transfected lysate using GPN1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with XAB1 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619