

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GSTA1-228


Anti GSTA1

Product type	Primary antibodies
Immunogen	GSTA1 (NP_665683.1, 1 a.a. ~ 222 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAEKPKLHYFNARGRMESTRWLLAAAGVEFEKFIKSAEDLDKLRNDGYLMFQQVPMVEI DGMKLVQTRAILNYIASKYNLYGKDIKERALIDMYIEGIADLGEMILLPVCPPPEEKDAK LALIKEKIKNRYFFPAFEKVLKSHGQDYLVGNKLSRADIHLVELLYYVEELDSSLISFPL LKALKTRISNLPTVKKFLQPGSPRKPPMDEKSLEEARKIFRF

Gene Alias cST2, GSTA1-1, GTH1, MGC131939

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~25.60 kDa using GSTA1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of GSTA1 expression in transfected 293T cell line by GSTA1 rabbit polyclonal antibody.


Lane 1: GSTA1 transfected lysate(25.60 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of GSTA1 transfected lysate

Immunoprecipitation



Immunoprecipitation of GSTA1 transfected lysate using GSTA1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GSTA1 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619