

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GTF2F2-217


Anti GTF2F2

Product type	Primary antibodies
Immunogen	GTF2F2 (AAH01771.1, 1 a.a. ~ 250 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAERGELDTGAKQNTGVWLVKVPKYLSQQWAKASGRGEVGLRIAKTQGRTEVSFTLNE DLANIHDIGGKPPASVSAPREHPFVLQSVGGQTLTVFTESSSDKLSLEGIVVQRAECRPAA SENYMRLKRLQIEESSKPVRLSQQLDKVTTNYKPVANHQYNIEYERKKKEGGKRARADK QHVLDMLFSAFEKHQYYNLKDLVDITKQPVVYLKEILKEIGVQNVKGIHKNTWELKPEYR HYQGEEKSD*

Gene Alias	ATP2B2, DKFZp686G08106, DKFZp686M088, MXRA1, PMCA4, PMCA4b, PMCA4x
-------------------	--

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~129.40 kDa using ATP2B4 transfected 293T lysate..
--------------------------	--

Western Blotting (Transfected lysate)



Western Blot analysis of GTF2F2 expression in transfected 293T cell line by GTF2F2 rabbit polyclonal antibody.

Lane 1: GTF2F2 transfected lysate(27.50 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619