

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GTF2H4-361

Anti GTF2H4

Product type	Primary antibodies
Immunogen	GTF2H4 (NP_001508.1, 1 a.a. ~ 462 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MESTPSRGLNRVHLQCRNLQEFLGGLSPGVLDRLYGHPATCLAVFRELP SLAKNWVMRML FLEQPLPQAAVALWVKKEFSKAQEESTGLLSGLRIWHTQLLPGGLQGLILNPIFRQNLRI ALLGGGKAWSDDTSQLGPDKHARDVPSLDKYAEERWEVVLHFMVGSPSAAVSQDLAQLLS QAGLMKSTEPGEPPCITSAGFQFLLLDTPAQLWYFMLQYLQTAQSRGMDLVEILSFLFQL SFSTLGKDYSVEGMSDSSLNFLQHLREFGLVFQKRKSRYYPTRLAINLSSGVSGAGGT VHQPFGFIVVETNYRLYAYTESELQIALIALFSEMLYRFPNMVVAQVTRESVQQAIASGIT AQQIIHFLRTRAHPVMLKQTPVLPPTITDQIRLWELERDRLRFTEGVLYNQFLSQVDFEL LLAHARELGVLVFENSAKRLMVVTPAGHSDVKRFWKRQKHSS


TFB2, TFIH

Application notes

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of GTF2H4 transfected lysate.

Immunoprecipitation



Immunoprecipitation of GTF2H4 transfected lysate using GTF2H4 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GTF2H4 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619