

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GUK1-292

Anti GUK1

Product type Primary antibodies

Immunogen GUK1 (NP_000849.1, 1 a.a. ~ 197 a.a) full-length human protein..

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit


Immunogen Sequence MSGPRPVVLSGPGAGKSTLLKRLQEHSGIFGFSVSHTRNPRPGEENGKDYYFVTREV
MQRDIAAGDFIEHAEFSGNLYGTSKVAVQAVQAMNRCVLDVDLQGVNRNIKATDLRPIYI
SVQPPSLHVLEQRLRQRNTETEEESLVKRLAAQAQADMESSEKEPGLFDVVIINDSLDQAYAE
LKEALSEEIKKAQRTGA

Gene Alias GMK

Application WB

notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~21.70 kDa using GUK1 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of GUK1 expression in transfected 293T cell line by GUK1 rabbit polyclonal antibody.

Lane 1: GUK1 transfected lysate(21.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of GUK1 expression in human spleen (Predicted Mw: 21.5940 kDa).

Western Blotting (Tissue lysate)




GUK1 rabbit polyclonal antibody. Western Blot analysis of GUK1 expression in human spleen.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of GUK1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of GUK1 transfected lysate using GUK1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GUK1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619