

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-GUSB-233

Anti GUSB

Product type	Primary antibodies
Immunogen	GUSB (NP_000172.1, 1 a.a. ~ 651 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MARGSAVAWAALGPLLWGCALGLQGGMLYPQESPSRECKELDGLWSFRADFSNRRRGFE EQWYRRPLWESGPTVDMVPVPSSFNDISQDWRLRHFGVWVWYEREVILPERWTQDLRTRVV LRIGSAHSYAIVWVNGVDLEHEGGYLPFEADISNLVQVGPLPSRLRITIAINNTLTPTT LPPGTIQYLTDTSKYPKGYFVQNTYFDFFNYAGLQRSVLLYTTPTTYIDDITVTTTSVEQD SGLVNYQISVKGSNLFKLEVRLLDAENKVAVANGTGTQGQLKVPGVSLWWPYLMHERPAYL YSLEVQLTAQTSGLPVSDFYTLPGVIRTVAVTKSQFLINGKPFYFHGVNKHEDADIRGKG FDWPLLVKDFNLLRWLGANAFRTSHYPYAEVVMQMCDRYGIVVIDECPGVGLALPQFFNN VSLHHHMQVMEEVVRDKNHPAVVMWSVANEPASHLESAGYYLKMVIAHTKSLDPSRPVT FVSNSNYAADKGAPYVDVICLSYYSWYHDYGHLELIQLQLATQFENWYKKYQKPIIQSE YGAETIAGFHQDPPLMFTEEYQKSLLEQYHLGLDQKRRKYVVGELIWNFADFMTESQSPTR VLGNKKGIFTRQRQPKSAAFLLRERYWKIANETRYPHSVAKSQCLENSPFT
Gene Alias	FLJ39445, MPS7
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~74.70 kDa using GUSB transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of GUSB expression in transfected 293T cell line by GUSB rabbit polyclonal antibody.

Lane 1: GUSB transfected lysate(74.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619