

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-HAO1-282

Anti HAO1

Product type Primary antibodies

Immunogen HAO1 (NP_060015.1, 1 a.a. ~ 370 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MLPRLICINDYEQHAKSVLPKSIYDYYRSGANDEETLADNIAAFSRWKLYPRMLRNVAET
DLSTSVLGQRVSMPICVGATAMQORMAHVDGELATVRACQSLGTGMMMLSSWATSSIEEVAE
AGPEALRWLQLYIYKDREVTKKLVRQAEEKMGYKAIFVTVDTPYLGNRLLDDVRNRFKLPPQ
LRMKNFETSTLSFSPEENFGDDSGLAAYVAKAIDPSISWEDIKWLRRRLTSLPIVAKGILR
GDDAREAVKHGLNGILVSNHGARQLDGVPATIDVLPEIVEAVEGKVEVFLDGGVVRKGTDV
LKALALGAKAVFVGRPIVWGLAFQGEKGVQDVLEILKEEFRLAMALSGCQNVKVIDKTLV
RKNPLAVSKI

Gene Alias GOX, GOX1, HAOX1, MGC142225, MGC142227

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.90 kDa using HAO1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of HAO1 expression in transfected 293T cell line by HAO1 rabbit polyclonal antibody.

Lane 1: HAO1 transfected lysate(40.90 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619