

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-HLCS-242

Anti HLCS

Product type Primary antibodies

Immunogen HLCS (NP_000402.3, 1 a.a. ~ 726 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit


Immunogen Sequence MEDRLHMDNGLVPQKIVSVHLQDSTLKEVKDQVSNKQAQILEPKPEPSLEIKPEQDGMEH
VGRDDPKALGEEKQRRGSASGSEPAGDSDRGGGPVEHYHLHLSSCHECLELENSTIESV
KFASAENIPDLPHYDYSSSLESVADETSPEREGRVNLTGKAPNILLYVGSDSQEALGRFH
EVRSLADCVDIDSYILYHLLDSALRDPWTDNCLLLVIATRESIPEDLYQKFMAYLSQG
GKVLGLSSSFTFGGFQVTSKGALHKTQVQNLVFSKADQSEVKLSVLSSGCRYQEGPVRLSP
GRLQGHLENEDKDRMIVHVPFGTRGGEAVLCQVHLELPPSSNIVQTPEDFNLLKSSNFRR
YEVLRILTTLGLSCDMKQVPALTPLYLLSAAEEIRDPLMQWLKGKVDSEGEIKSGQLSL
RFVSSYVSEVEITPSCIPVVTNMEAFSSEHFNLEIYRQNLQTKQLGKVLFAEVTPTTMR
LLDGLMFQTPQEMGLIVIAARQTEGKGRGGNVWLSVPGCALSTLLISIPRSQLGQRIPF
VQHLMSVAVVEAVRSIPEYQDINLRVKWPNDIYYSDLMKIGGVLVNSTLMGETFYILIGC
GFNVTNSNPTICINDLITEYNKQHKAEKPLRADYLIARVVTVLEKLIKEFQDKGPNSVL
PLYRYRWVHSGQQVHLGSAEGPKVSIVGLDDSGFLQVHQEGGEVTVHPDGNSFDMRLNL
ILPKRR

Gene Alias HCS

Application WB

notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~80.80 kDa using HLCS transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of HLCS expression in transfected 293T cell line by HLCS rabbit polyclonal antibody.


Lane 1: HLCS transfected lysate(80.80 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of HLCS transfected lysate.

Immunoprecipitation



Immunoprecipitation of HLCS transfected lysate using HLCS rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with HLCS rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619