



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-HMGCR-161

Anti HMGCR

Product type	Primary antibodies
Immunogen	HMGCR (AAH33692.1, 1 a.a. ~ 835 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit


Immunogen Sequence

MLSRLFRMHGLFVASHPWEVIVGTVTLTICMMSMMFTGNNKICGWNYECPKFEEDVLSS
 DIIILTITRCIAILYIFQFQNLRLQSGKYILGIAGLFTIFSSFVFSTVVIHFLDKELTG
 LNEALPFFLLIDLSRASTLAKFALSSNSQDEVRENIARGMAILGPTFTLDALVECLVIG
 VGTMSGVRQLEIMCCFGCMSVLANYFVFMTFFPACVSLVLELSRESREGRPIWQLSHFAR
 VLEEEENKPNPVTQRVKMIMSLGLVLVHAHSRWIADPSPQNSTADTSKVSLGLDENVSKR
 IEPVSLWQFYLSKMISMDIEQVITLSLALLAVKYIFFEQTETESTLSLKNPITSPVVT
 QKKVPDNCRRREPMLVRNMQKCDSEVEETGINRERKVEVIKPLVAETDTPNRATFVVGNS
 SLLDTSSVLVTQEPEIELPREPRNEECLQILGNAEKGAKFLSDAEIQLVNAKHIPAYK
 LETLMETHERGVSIRRQLLSKKLSEPSLQYLPYRDYNYSLGGGASSRVLADGMTRGPV
 VRLPRACDSEVKAWLETSEGFAVIKEAFDSTSRFARLQKLHTSIAGRNLIRFQSRSGD
 AMGMNMISKGTEKALS KLHEYFPEMQILAVSGNYCTDKKPAAINWIEGRGKSVVCEAVIP
 AKVVREVLKTTTEAMIEVNINKNLVGSAMAGSIGGYNAHAANIVTAIYIACGQDAAQNVG
 SSNCITLMEASGPTNEDLYISCTMPSIEIGTVGGGTNLLPQQAQLQMLGVQGACKDNPGE
 NARQLARIVCGTVMAGELSLMAALAAGHLVKSHMIHNRSKINLQDLQGACTKKTA

Gene Alias

Application notes **IP**
 Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of HMGCR transfected lysate

Immunoprecipitation



Immunoprecipitation of HMGCR transfected lysate using HMGCR rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with HMGCR purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619