



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-IL22-130

Anti IL-22

Product type	Primary antibodies
Immunogen	IL22 (NP_065386.1, 1 a.a. ~ 179 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAALQKSVSSFLMGTLATSCLLLLALLVQGGAAAPISSHCRLDKSNFQQPYITNRTFMLA KEASLADNNTDVRLLIGEKLFGVSMSERCYLMKQVLNFTLEEVLPQSDRFQPYMQEVVP FLARLSNRLSTCHIEGDDLHIQRNVQKLKDTVKKLGESGEIKAIGELDLLFMSLRNACI
Gene Alias	IL-21, IL-22, IL-D110, IL-TIF, IL21, ILTIF, MGC79382, MGC79384, TIFIL-23, TIFa, zcyto18
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~20.00 kDa using IL22 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of IL22 expression in transfected 293T cell line by IL22 rabbit polyclonal antibody.

Lane 1: IL22 transfected lysate (20.00 kDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619