



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-I P-151

Anti IL10RB


Product type	Primary antibodies
Immunogen	IL10RB (AAH01903.1, 1 a.a. ~ 325 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MAWSLGSWLGGCLLVSA LGMVPPPENVRMNSVNFKNILQWESPAFAKGNLTFTAQYLSYR IFQDKCMNTTLTECDFSSLSKYGDHTLRVRAEFADEHSDWVNITFCPVDDTIIGPPGMQV EVLADSLHMRFLAPKIENEYETWTMKNVYNSWTYNVQYWKNGTDEKFQITPQYDFEVLRN LEPWTTYCVQVRGFLPDRNKAGEWSEPVCEQTTHDETVP SWMVAVILMASVFMVCLALLG CFALLWCVYKKT KYAFSPRNSLPQHLKEFLGHPHHNTLLFFSFPLSDENDVFDKLSVIAE DSESGKQNP GDSCSLGT PPGQPQS
--------------------	---

Gene Alias	CDW210B, CRF2-4, CRFB4, D21S58, D21S66, IL-10R2
------------	---

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.75 kDa using IL10RB transfected 293T lysate
-------------------	---

Western Blotting (Transfected lysate)



Western Blot analysis of IL10RB expression in transfected 293T cell line by IL10RB rabbit polyclonal antibody.


Lane 1: IL10RB transfected lysate(35.75 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of IL10RB transfected lysate.

Immunoprecipitation



Immunoprecipitation of IL10RB transfected lysate using IL10RB rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with IL10RB mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619