

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- IP-171

Anti NAPRT1

Product type	Primary antibodies
Immunogen	NAPRT1 (NP_660202.2, 1 a.a. ~ 466 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MALGYWRAGRARDAAEFELFFRRCPFGGAFALAAGLRDCVRFLRAFRLRDADVQFLASVL PPDTPDAFFEHLRALDCSEVTVRALPEGSLAFPGVPLLQVSGPLLVVQLLETPLLCLVSY ASLVATNAARLRLIAGPEKRLLEMGLRRAQGPDGGLTASTYSYLGGFSSSNVLAGQLRG VPVAGTLAHSFVTSFSGSEVPPDPMLAPAAGEGPGVDLAAKAQVWLEQVCAHLGLGVQEP HPGERAAFVAYALAFPRAFQGLLDITYSVWRSGLPNFLAVALALGELGYRAVGVRDLSGDL LQQAQEIIRKVFRAAAQFQVPWLESVLIVVSNNIDEEALARLAQEGSEVNVIGIGTSVVT CPQQPSLGGVYKLVAVGGQPRMKLTEDPEKQTLPGSKAAFRLLGSDGSPLMDMLQLAEEP VPQAGQELRVWPPGAQEPCTVRPAQVEPLLRLCLQQGQVAAPPPSS
Gene Alias	PP3856

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.60 kDa using NAPRT1 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of NAPRT1 expression in transfected 293T cell line by NAPRT1 rabbit polyclonal antibody.


Lane 1: NAPRT1 transfected lysate(49.60 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of NAPRT1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of NAPRT1 transfected lysate using NAPRT1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with NAPRT1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619