



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-180


Anti LGALS3BP

Product type	Primary antibodies
Immunogen	LGALS3BP (NP_005558.1, 1 a.a. ~ 585 a.a) full-length human protein
Form	Liquid
Volume	100 uL
Specificity	Human
Host	Rabbit
Immunogen Sequence	MTPPRLFWVWLLVAGTQGVNDGDMRLADGGATNQGRVEIFYRGQWGTVCNLDLTDASV VCRALGFENATQALGRAAFGQGSGPIMLDEVQCTGTEASLADCKSLGWLKSNCRHERDAG VVCTNETRSTHTLDLSRELSEALGQIFDSQRGCDLSISVNVQGEDALGFCGHTVILTANL EAQALWKEPGSNVTMSVDAECVPMVRDLLRYFYSSRRIDITLSSVKCFHKLASAYGARQLQ GYCASLFAILLPQDPSFQMPLDLYAYAVATGDALLEKLCLQFLAWNFEALTQAEAWPSVP TDLLQLLLPRSDLAVPSELALLKAVDTWSWGERASHEEVEGLVEKIRFPMMLPEELFELQ FNLSLYWSHEALFQKKTLLQALEFHTVPFQLLARYKGLNLTEDTYKPRIYTSPTWSAFVTD SSWSARKSQLVYQSRRGPLVKYSSDYFQAPSDYRYYPYQSFQTPQHPSFLFQDKRVSWSL VYLPTIQSCWNYGFSCSSDELPVLGLTKSGGSDRTIAYENKALMLCEGLFVADVTDFEGW KAAIPSAALDTNSSKSTSSFPAGHFNGFRTVIRPFYLTNSSGVD

Gene Alias 90K, BTBD17B, MAC-2-BP

Application IP
notes Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LGALS3BP transfected lysate

Immunoprecipitation



Immunoprecipitation of LGALS3BP transfected lysate using LGALS3BP rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LGALS3BP purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610
FAX : +81-3-5632-9618 F AX : +81-3-5632-9619