



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-196

Anti PHF10

Product type	Primary antibodies
Immunogen	PHF10 (NP_060758.1, 1 a.a. ~ 410 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MLQEQVSEYLGVTSEFKRKYPDLERRDLSHKEKLYLRELVITETQCTLGLTALRSDEVID LMIKEYPAKHAEYSVILQEKERQRITDHYKEYSQMQQQNTQKVEASKVPEYIKKAAKAA EFNSNLNRERMEERRAYFDLQTHVIQVPQGKYKVLPTERTKVSSYPVALIPGQFQEYYKR YSPDELRYPPLNTALYEPPLDPELPALDSDGSDDGEDGRGDEKRKNKGTSDSSSGNVSE GESPPDSQEDSFQGRQKSKDKAATPRKDGPKRSVLSKSVPGYKPKVIPNAICGICLKGE SNKKGKAESLIHCSQCENSGHPSCLDMTMELVSMIKTYPWQCMECKTCICGQPHHEEEM MFCDMCDRGYHTFCVGLGAIPSGRWICDCCQRAPPTPRKVGRRGKNSKEG
Gene Alias	FLJ10709
Application notes	IP Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PHF10 transfected lysate

Immunoprecipitation



Immunoprecipitation of PHF10 transfected lysate using PHF10 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PHF10 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619