

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-224

Anti SULT2B1

Product type	Primary antibodies
Immunogen	SULT2B1 (NP_814444.1, 1 a.a. ~ 365 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MDGPAEPQIPGLWDTYEDDISEISQKLPGEYFRYKGVFPFVGLYSLESISLAENTQDV RD DDIFIITYPKSGTTWMIEIICLILKEGDPWIRSVPIWERAPWCETIVGAFSLPDQYSPR LMSSHLPIQIFTKAFFSSKAKVIYMGRNPRDVVSLYHYSKIAGQLKDPGTPDQFLRDFL KGEVQFGSWFDHIKGWLRMKGKDNFLFITYEELQQDLQGSVERICGFLGRPLGKEALGSV VAHSTFSAMKANTMSNYTLLPPSLLDHRGAFLRKGVCGDWKNHFTVAQSEAFDRAYRKQ MRGMPTFPWDEDPEEDGSPDPEPSPEPEPKPSLEPNTSLEREPRPNSSPSPSPGQASETP HPRPS
Gene Alias	HSST2
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~41.30 kDa using SULT2B1 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of SULT2B1 expression in transfected 293T cell line by SULT2B1 rabbit polyclonal antibody


Lane 1: SULT2B1 transfected lysate(41.30 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of SULT2B1 transfected lysate

Immunoprecipitation



Immunoprecipitation of SULT2B1 transfected lysate using SULT2B1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with SULT2B1 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619