

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- IP-279

Anti AADAT

Product type	Primary antibodies
Immunogen	AADAT (AAH31068.1, 1 a.a. ~ 429 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MNYARFITAASAARNPSPIRTMSEKRADILSRGPKSMISLAGGLPNPNMFPPKTAVITVE NGKTIQFGEEMMKRALQYSPSAGIPELLSWLKQLQIKLHNPPTIHYQPSQGQMDLCVTSG SQQGLCKVFEMIINPGDNLVLLDEPAYSGTLQSLHPLGCNIINVASDESGIVPDSLRLDILS RWKPEDAKNPQKNTPKFLYTPNGNPTGNSLTSEKKEIYELARKYDFLIIEDDPYYFL QFNKFRVPTFLSMDVDGRVIRADSFSKIISGLRIGFLTGPKPLIERVILHIQVSTLHPS TFNQLMISQLLHEWGEEGFMAHVDRVIDFYSNQKDAILAAADKWLTGLAEWHVPAAGMFL WIKVKGINDVKELIEEKAVKMGVLMPLPGNAFYVDSSAPSPYLRSFSSASPEQMDVAFQV LAQLIKESL

Gene Alias	KAT2, KATII
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~47.90 kDa using AADAT transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of AADAT expression in transfected 293T cell line by AADAT rabbit polyclonal antibody.


Lane 1: AADAT transfected lysate(47.90 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of AADAT transfected lysate.

Immunoprecipitation



Immunoprecipitation of AADAT transfected lysate using AADAT rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with AADAT purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619