

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-362

Anti PFKFB2

Product type Primary antibodies

Immunogen PFKFB2 (NP_006203.2, 1 a.a. ~ 505 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MSGASSSEQNNNSYETKTPNLRMSEKKCSWASYMTNSPTLIVMIGLPARGKTYVSKKLTR
YLNWIGVPTKVFNLGVYRREAVKSYKSYDFFRHDNEEAMKIRKQCALVALEDVKAYLTEE
NGQIAVFDATNTTTRERRDMILNFAEQNSFKVFFVESVCDDPDVIAANILEVKVSSPDYPE
RNRENVMEDFLKRIECYKVTYRPLDPDNYDKDLSFIKVINVGQRFLVNRVQDYIQSKIVY
YLMNIHVQPRTIYLCRHGESEFNLLGKIGGDSGLSVRGKQFAQALRKFLLEEQEITDLKVV
TSQLKRTIQTAEISLGVPEQWKILNEIDAGVCEEMTYAEIEKRYPEEFALRDQEKYLYRY
PGGESYQDLVQRLEPVMELERQGNVLVISHQAVMRCLLAYFLDKGADELPYLRCPLHTI
FKLTPVAYGCKVETIKLNVEAVNTHRDKPTNFPKNQTPVRMRRNSFTPLSSSNTIRRP
NYSVGSRLPLKPLSPLRAQDMQEGAD


Gene Alias DKFZp781D2217, MGC138308, MGC138310, PFK-2/FBPase-2

Application notes

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PFKFB2 transfected lysate

Immunoprecipitation



Immunoprecipitation of PFKFB2 transfected lysate using PFKFB2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PFKFB2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619