

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-371

Anti SOHLH2

Product Primary antibodies

type

Immunogen SOHLH2 (AAH25383.1, 1 a.a. ~ 225 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MASSIICQEHQCISGQAKIDILLVGDVTVGYLADTVQKLFANIAEVTITISDTKEAAALL
DDCIFNMVLLKVPSSLSAEELEAIKLIRFGKKKNTHSLFVFIIPENFKGCISGHGMDIAL
TEPLTMEKMSNVVKYWTTCPSTNTVKTENATGPEELGLPLQRSYSEHLGYFPTDLFACSES
LRNGNGLELNASLSEFEKNKKISLLHSSKEKLRRRLYRKHSSFCFW


FLJ20449, TEB1, bA121N13.2

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~25.10 kDa using SOHLH2 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of SOHLH2 expression in transfected 293T cell line by SOHLH2 rabbit polyclonal antibody.

Lane 1: SOHLH2 transfected lysate(25.10 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of SOHLH2 expression in human liver
(Predicted Mw: 46.9410 kDa).

Western Blotting (Tissue lysate)




SOHLH2 rabbit polyclonal antibody. Western Blot analysis of SOHLH2 expression in human liver.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of SOHLH2 transfected lysate

Immunoprecipitation



Immunoprecipitation of SOHLH2 transfected lysate using SOHLH2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TEB1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer.
Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619