

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-477


Anti CRYZ

Product type	Primary antibodies
Immunogen	CRYZ (NP_001880.2, 1 a.a. ~ 329 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MATGQKLMRAVRVFEFGGPEVLKLRSDIAVPIPKDHQVLIKVHACGVNPVETYIRSGTYS RKPLLPTYTPGSDVAGVIEAVGDNASAFKKGDRVFTSSTISGGYAEYALAADHTVYKLPEK LDFKQGAAGIPYFTAYRALIHSACVKAGESVLVHGASGGVGLAACQIARAYGLKILGTA GTEEGQKIVLQNGAHEVFNHREVNYIDKIKKYVGEKGDIIIEMLANVNL SKDLSLLSHG GRVIVVGSRGTEINPRDTMAKESSIIGVTLFSSTKEEFQQYAAALQAGMEIGWLKPVIG SQYPLEKVAEAHENIIHGSGATGKMILL

Gene Alias DKFZp779E0834, FLJ41475

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.20 kDa using CRYZ transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of CRYZ expression in transfected 293T cell line by CRYZ rabbit polyclonal antibody.


Lane 1: CRYZ transfected lysate(35.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of CRYZ transfected lysate

Immunoprecipitation



Immunoprecipitation of CRYZ transfected lysate using CRYZ rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with CRYZ purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619