

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-491

Anti HMGB1

Product type Primary antibodies

Immunogen HMGB1 (NP_002119.1, 1 a.a. ~ 215 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MGKGDPPKPRGKMSSYAFFVQTCREEHKKKHDPASVNFSEFSKKCSERWKTMSAKEKGKF
EDMAKADKARYEREMKTYIPPKGETKKKFKDPNAPKRPPSAFFLCSEYRPKIKGEHPGL
SIGDVAKKLGEMWNNNTAADDKQPYEKKAALKKEYEKDIAAYRAKGKPDAAKKGVVKA EK
SKKKKEEEDEEDEDEDEDEDEDEDEDEDEDEDEDEDEDEDDDE

Gene Alias DKFZp686A04236, HMG1, HMG3, SBP-1

Application notes **IP**
Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of HMGB1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of HMGB1 transfected lysate using HMGB1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with HMGB1 monoclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619