

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-496

Anti NAB2

Product type	Primary antibodies
Immunogen	NAB2 (AAH07756.1, 1 a.a. ~ 241 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MHRAPSPTAEQPPGGGDSARRTLQPRLKPSARAMALPRTLGEQLYRVLQRANLLSYYET FIQQGGDDVQQLCEAGEEEFLEIMALVGMATKPLHVRRLQKALREWATNPGLFSQPVP PVSSIPLFKISETAGTRKGSMSNGHSGPGEKAGSARSFSPKSPLELGEKLSPLPGGPGAG DPRIWPGRSTPESDVGAGGEEEEAGSRWDWGWRSRPTGARDGTPACLGRVWMDICRLWGHVQ G


Gene Alias MADER, MGC75085

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~25.90 kDa using NAB2 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of NAB2 expression in transfected 293T cell line by NAB2 rabbit polyclonal antibody.


Lane 1: NAB2 transfected lysate(25.90 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of NAB2 transfected lysate

Immunoprecipitation



Immunoprecipitation of NAB2 transfected lysate using NAB2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with NAB2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619