

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-538

Anti FARS2

Product type Primary antibodies

Immunogen FARS2 (NP_006558.1, 1 a.a. ~ 451 a.a) full-length human protein..

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MVGSALRRGAHAYVYLVS KASHISRGHQHQAWGSRPPAAECATQRAPGSVVELLGKSY PQ
DDHSNLTRKVLTRVGRNLHNQQHHPLWLIKERVKEHFYKQYVGRFGTPLFSVYDNLSPV
TTWQNFDSLLIPADHPSRKKGDNYLNRTHMLRAHTSAHQWDL LHAGLDAFLVVG DVYRR
DQIDSQHYPIFHQLEAVRLF SKHEL FAGIKDGESLQLFEQSSRSAHKQETH TMEAVKLVE
FDLKQTLTRLMAHLFGDELEIRWVDCYFPFTHPSFEMEINFHGEWLEVLGCGVMEQQLVN
SAGAQDRIGWAFGLGLERLAMILYDIPDIRLFWCEDERFLKQFCVSNINQKV KFQPLSKY
PAVINDISFWLPSENYAENDFYDLVRTIGGDLVEKVDLIDKFVHPKTHKTS HCYRITYRH
MERTLSQREVRHIHQALQEAAVQLLGVEGRF


Gene Alias FARS1, HSPC320, PheRS, dJ520B18.2

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~52.40 kDa using FARS2 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of FARS2 expression in transfected 293T cell line by FARS2 rabbit polyclonal antibody.


Lane 1: FARS2 transfected lysate(52.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of FARS2 transfected lysate

Immunoprecipitation



Immunoprecipitation of FARS2 transfected lysate using FARS2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with FARS2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619