

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-556

Anti CDCA8


Product type	Primary antibodies
Immunogen	CDCA8 (NP_060571.1, 1 a.a. ~ 280 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAPRKGSSRVAKTNSLRRRKLASFLKDFDREVEIRIKQIESDRQNLLKEVDNLYNIEILR LPKALREMNWLDYFALGGNKQALEEAATADLDITEINKLTAEAIQTPLKSAKTRKVIQVD EMIVEEEEEENERKNLQTARVKRCPPSKKRTQSIQGKGKGKRSSRANTVTPAVGRLEVS MVKPTPGLTPRFDSRVFKTPGLRTPAAGERIYNISGNGSPLADSKEIFLTVPVGGGESLR LLASDLQRHSIAQLDPEALGNIKKLSNRLAQICSSIRTHK

Gene Alias BOR, BOREALIN, FLJ10468, FLJ12042

Application notes

WB
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~31.30 kDa using CDCA8 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CDCA8 expression in transfected 293T cell line by CDCA8 rabbit polyclonal antibody.


Lane 1: CDCA8 transfected lysate(31.30 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of CDCA8 transfected lysate

Immunoprecipitation



Immunoprecipitation of CDCA8 transfected lysate using CDCA8 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with CDCA8 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619