

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-557

Anti ASXL2

Product type Primary antibodies

Immunogen ASXL2 (AAH42999.1, 1 a.a. ~ 918 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MKRTKCADIDVETPDSILVNTNLRALINKHTFSVLPGDCQQRLLLLLPEVDRQVGPDGLM
KLNGSALNNEFFTSAAQGWKERLSEGEFTPEMQVRIRQEIEKEKKVEPWKEQFFESYYGQ
SSGLSLEDISKLTASPSDPKVKKTPAEQPKSMPVSEASLIRIVPVVSQSECKEEALQMSS
PGRKEECESQGEVQPNFSTSSEPLLSSALNTHLSSILPIKCPKDEDLLEQKPVTSAEQE
SEKNHLTTASNYNKSESQESLVTSPSKPKSPGVEKPIVKPTAGAGPQETNMKEPLATLVD
QSPESLKRKSSLTQEEAPVSWEKRPRTVENRQHQQPFQVSPQPFLNRGDRIQVRKVPPLK
IPVSRISPMFHPSPQVSPRARFPVSITSPNRTGARTLADIKAKAQLVKAQRAAAAAAAAAA
AAAASVGGTIPGPGPGGGGQPGEGGEGQTARGGSPGSDRVSETGKGPTLELAGTGSRRGT
RELLPCGPETQPQSETKTTPSQAQPHSVSGAQLQQTTPVPPTPAVSGACTSVPSPVHIEK
LDNEKLNPTTRATATVASVSHPPQGPSSCRQEKAQSPSTGFQLEDISTSQRFMLGFAGRRTSK
PAMAGHYLLNISTYGRGSEFRTHSVNPEDRFCLSSPTEALKMGYTDCKNATGESSSSK
EDDTDEESTGDEQESVTVKEEPQVSQSAGKGDTS SGPHSRETSTSDCLASKNVKAEIPL
NEQTTLSKENYLFTRGQTFDEKTLARDLIQAAQKQMAHAVRGKAIRSSPELFSSTVLPLP
ADSPTHQPLLLPPLQTPKLYGSPTQIGPSYRGMINVSTSSDMDHNSAVPGSQVSSNVGDV
MSFSVTVTIPASQAMNPSSHGQTIPVQAFSEENSIEGTPSKCYCRLKAMIMCKGCGAFC
HDDCIGPSKLCVSLVVR


Gene Alias ASXH2, DKFZp686C1968, FLJ10898, KIAA1685

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~98.40 kDa using ASXL2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ASXL2 expression in transfected 293T cell line by ASXL2 rabbit polyclonal antibody.


Lane 1: ASXL2 transfected lysate(98.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of ASXL2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of ASXL2 transfected lysate using ASXL2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with ASXL2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619