

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-KIAA1199-208

Anti KIAA1199

Product type	Primary antibodies
Immunogen	KIAA1199 (AAH20256.1, 1 a.a. ~ 992 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MGAAGRQDFLFKAMLTISWLTLTCTPGATSTVAAGCPDQSPQLPWNPGHDQDHHVHIGQ GKTLLLTSSATVYSIHISEGGKLVKDHDEPIVLRTRHILIDNGGELHAGSALCPFQGNF TIILYGRADEGIQPDPIYGLKYIGVGKGALELHGQKKLSWTFNLKTLHPGGMAEGGYFF ERSWGHHRGVIVHVIDPKSGTVIHSDRFDYRSKKESERLVQYLNAVDPGRILSVAVNDEG SRNLDDMARKAMTKLGSKHFLHLGFRHPWSFLTVMKGNPSSSVEDHIEYHGHRGSAAARVF KLFQTEHGEYFNVSLSSSEWVQDVEWTEWFDHDKVSQTKGGEKISDLWKAHPGKICNRPID IQATTMDGVNLSTEVVYKKGQDYRFACYDRGRACRSYRVFLCGKPVPRKLTVTIDTNVN STILNLEDNVQSWKPGDTLVIASTDYSMYQAEFQVLPCRSCAPNQVKVAGKPMYHIGE EIDGVDMAEVLGSLRNIIVMGEMEDKCPYRNHICNFFDFDTFGGHIKFAFGKAAHLE GTELKHMGGQLVGQYPIHFHLAGDVDERGGYDPPTYIRDLSIHHTFSRCVTVHGSNGLLI KDVVGYNSLGHCFFTEDGPEERNTFDHCLGLLVKSGTLLPSDRDSKMCKMITEDSYPGYI PKPRQDCNAVSTFWMANPNNNLINCAAAGSEETGFWFIFHHVPTGPSVGMYSPPGYSEHIP LGKFYNNRAHSNYRAGMIIDNGVKTTEASAKDKRPFLSIISARYSPHQDADPLKPREPAI IRHFIAYKNQDHGAWLRGGDVWLDSCRFAENGIGLTLASGGTFPYDDGSKQEIKNSLFVG ESGNVGTEMMDNRIWGPGLDHSGRITLPIGQNFPIRGILYDGPINIQNCTFRKFVALEG RHTSALAFRLNNAWQSCPHNNVTGIAFEDVPITSRVFFGEPGPWFNQLDMDGDKTSVFHD VDGSVSEYPGSYLTKNCKWHSLSKASGPGSG


Gene Alias TMEM2L

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~110.40 kDa using KIAA1199 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of KIAA1199 expression in transfected 293T cell line by KIAA1199 rabbit polyclonal antibody.


Lane 1: KIAA1199 transfected lysate(110.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of KIAA1199 transfected lysate.

Immunoprecipitation



Immunoprecipitation of KIAA1199 transfected lysate using KIAA1199 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with KIAA1199 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619