

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-KRT12-120

Anti KRT12

Product type	Primary antibodies
Immunogen	KRT12 (AAI56642.1, 1 a.a. ~ 494 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MDLSNNTMSLSVRTPGLSRRLSSQSVIGRPRGMSASSVSGSYGGSAFGFGASC GG GFSAASMF GSSSGFGGGSGSSMAGGLGAGYGRALGGGSFGGLGMGFGGSPGGGSLGILSGNDGGLSGSEKETMQNLNDRLASYLDKVRAL EEA NTELENKIREWYETRGTGTADASQSDYSKYYPLIEDLRNKIISASIGNAQLLLQIDNARLAAEDFRMKYENELALRQGV EADINGLR RVLD ELTLTRTDLEMQIESLNEELAYMKKNHEDELQSFRVGGPGEVSVEMDAAPGVDLTRLLNDMRAQYETIAEQNRKD AEAWFIEKSGELRKEISTNTEQLQSSKSEVTDLRRAFQNL EIELQSQLAMKKSLED SLAEAGDYCAQLSQVQQLISNLEAQLLQVRADAERQNV DHQRLLNVKARLELEIETYRRLLDGEAQGDGLEESLFVTD SKSQAQSTDSSKDPTKTRKIKTVVQEMVNGEVVSSQVQEIEELM
Gene Alias	K12
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~54.34 kDa using KRT12 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of KRT12 expression in transfected 293T cell line by KRT12 rabbit polyclonal antibody.

Lane 1: KRT12 transfected lysate(54.34 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619