

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-LARP6-188

Anti LARP6


Product type	Primary antibodies
Immunogen	LARP6 (NP_060827.2, 1 a.a. ~ 491 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAQSGGEARPGPKTAVQIRVAIQEAEDVDELEDEEEGAETRGAGDPARYLSPGWGSASEE EPSRGHSGTTASGGENEREDLEQEWKPPDEELIKKLVDQIEFYFSDENLEKDAFLKHKVR RNKLGYSVSKLLTSFKKVKHLTRDWRTTAHALKYSVLELNEDHRKVRRTTPVPLFPNEN LPSKMLLVYDLYLSPKLWALATPQKNGRVQEKVMEHLLKLFGTFGVISSVRILKPGRELP PDIRRISSRYSQVGTQECAIVEFEEVEAAIKAHEFMITESQGKENMKAVLIGMKPPKKKP AKDKNHDEEPTASIHNLNKSINRVEELQYMGDESSANSSSDPESNPTSPMAGRRHAATNK LSPSGHQNLFLSPNASPCTSPWSSPLAQRKGVSRKSPLAEEGRLNCSTSPEIFRKCMDYS SDSSVTPSGSPWVRRRRQAEMGTQEKSPGTSPLLSRKMQTADGLPVGVLRLPRGPDNTRG FHGHERSRACV
Gene Alias	ACHN, FLJ11196

Application notes

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LARP6 transfected lysate

Immunoprecipitation



Immunoprecipitation of LARP6 transfected lysate using LARP6 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LARP6 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619