



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-LGALS3BP-180

Anti LGALS3BP


Product type	Primary antibodies
Immunogen	LGALS3BP (NP_005558.1, 1 a.a. ~ 585 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence
MTPPRLFWVWLLVAGTQGVNDGDMRLADGGATNQGRVEIFYRGQWGTVCDNLWDLTDASV
VCRALGFENATQALGRAAFGQSGPIMLDEVQCTGTEASLADCKSLGWLKSNCRHERDAG
VVCTNETRSTHTLDLSRELSEALGQIFDSQRGCDLSISVNVQGEDALGFCGHTVILTANL
EAQALWKEPGSNVTMSVDAECVPMVRDLLRYFYSSRRIDITLSSVKCFHKLASAYGARQLQ
GYCASLFAILLPQDPSFQMPLDLYAYAVATGDALLEKLCQFLAWNFEALTQAEAWPSVP
TDLLQLLPRSDLAVPSELALLKAVDTWSWGERASHEEVEGLVEKIRFPMMLPEELFELQ
FNLSLYWSHEALFQKKTLLQALEFHTVPFQLLARYKGLNLTEDTYKPRIYTSPTWSAFVTD
SSWSARKSQLVYQSRRGPLVKYSSDYFQAPSDYRYYPYQSFQTPQHPSFLFQDKRVSWSL
VYLPTIQSCWNYGFSCSSDELPLVGLTKSGGSDRTIAYENKALMLCEGLFVADVDFEGW
KAAIPALDTNSSKSTSSFPCPAGHFNGFRTVIRPFYLTNSSGVD

Gene Alias 90K, BTBD17B, MAC-2-BP

Application notes **IP**
Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LGALS3BP transfected lysate

Immunoprecipitation



Immunoprecipitation of LGALS3BP transfected lysate using LGALS3BP rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LGALS3BP purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619