

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-LRSAM1-139


Anti LRSAM1

Product type	Primary antibodies
Immunogen	LRSAM1 (NP_001005373.1, 1 a.a. ~ 723 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MPLFFRKRPSEEARKRLEYQMCLAKEAGADDILDISKCELSIEIPFGAFATCKVLQKKVL IVHTNHLTSLLPKSCSLLSLATIKVLDLHDNQLTALPDDLGLTALQVLNVERNQLMQLP RSIGNLTQLQTLNVKDNKLKELPDTVGE LRSRLTLNISGNEIQRLPQMLAHVRTLEMLS DASAMVYPPREVCAGAGTAAILQFLCKESGLEYYPPSQYLLPILEQDGIENS RDSPDGPTD RFSREELEWQNRFS DYEK RKEQKMLEKLEFERRLELGQREHTQLLQQSSSQKDEILQTVK EEQSRLEQGLSEHQ RHLNAERQRLQEQLKQTEQNISSRIQKLLQDNQRQKKSSEILKSLE NERIRMEQLMSITQEETESLRRRDVASAMQQMLTESCKNRLIQMAYESQRQNLVQQACSS MAEMDERFQQILSWQQMDQNK AISQILQESAMQKA AFEALQVKKDLMHRQIRSIKLIET ELLQLTQLELKRKSLDTESLQEMISEQRWALSSLLQQLLKEKQQREEELREILTELEAKS ETRQENYWLIIQYQRLLNQKPLSLKLQEEGMERQLVALLEELSAEHYLP IFAHHRLSLDLL SQMSPGDLAKVGVSEAGLQHEILRRVQELLDAARIQPELKPPMGEVVTPTAPQEPPESVR PSAPPAELEVQASECVVCLEREAQMIFLNCGHVCCCQCCQPLRTCPLCRQDIAQRLRIY HSS

Gene Alias

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~83.60 kDa using LRSAM1 transfected 293T lysate
--------------------------	---


Western Blotting (Cell lysate)



LRSAM1 rabbit polyclonal antibody. Western Blot analysis of LRSAM1 expression in HeLa.

Antibody dilution: 1:1000. Western Blot analysis of LRSAM1 expression in human liver (Predicted Mw: 83.5940 kDa).

Western Blotting (Tissue lysate)




LRSAM1 rabbit polyclonal antibody. Western Blot analysis of LRSAM1 expression in human liver.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LRSAM1 transfected lysate

Immunoprecipitation



Immunoprecipitation of LRSAM1 transfected lysate using LRSAM1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LRSAM1 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619