

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-MAN1A1-243**

# Anti MAN1A1

**Product type** Primary antibodies

**Immunogen** MAN1A1 (AAH65827.1, 1 a.a. ~ 367 a.a) full-length human protein.

**Form** Liquid  
Protein A purified

**Volume** 100ul

**Specificity** Human


**Host** Rabbit

**Immunogen Sequence** MPVGGLPLFSSPAGGVLGGGLGGGGGRKGS GPAALRLTEKFVLLLVFSAFITLCFGAIF  
FLPDSSKLLSGVLFHSSPALQPAADHKPGPGARAEDAAEGRARRREEGAPGDPEAALEDN  
LARIRENHERALREAKETLQKLPEEIQRDILLEKKKVAQDQLRDKAPFRGLPPVDFVPPI  
GVESREPADAAIREKRAKIKEEIEARKGQANCPGSSSSCEVEIQMMKHAWN NYKGYAWGL  
NELKPI SKGGHSSSLFGNIK GATIVDALDTLFIMEMKHEFEEAKSWVEENLDFNVNAEIS  
VFEVNIRFVGGLLSAYYLSGEEIFRKKAVELGVKLLPAFHTPSGIPWALLNMKRWMCNIP  
RAMEATI

**Gene Alias** HUMM3, HUMM9, MAN9

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.00 kDa using MAN1A1 transfected 293T lysate.

**Western Blotting (Transfected lysate)**




Western Blot analysis of MAN1A1 expression in transfected 293T cell line by MAN1A1 rabbit polyclonal antibody.

Lane 1: MAN1A1 transfected lysate(40.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of MAN1A1 expression in A-431 (Predicted Mw: 72.9690 kDa).

**Western Blotting (Cell lysate)**



MAN1A1 rabbit polyclonal antibody. Western Blot analysis of MAN1A1 expression in A-431.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619