

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-MAN1A2-239

Anti MAN1A2

Product type	Primary antibodies
Immunogen	MAN1A2 (NP_006690.1, 1 a.a. ~ 641 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MTTPALLPLSGRRIPPLNLGPPSFPHHRATLRLSEKFILLILSAFITLCFGAFFFLPDS SKHKRFDLGLDVLIPHVDAGKGAKNPGVFLIHGPDEHRHREEEERLRNKIRADHEKALE EAKEKLRSREEIRAEIQTEKNKVQEMKIKENKPLPPVIPNLVGIRGGDPEDNDIREK REKIKEMMKHAWDNYRTYGWHNELRPIARKGHSPNIFGSSQMGATIVDALDTLYIMGLH DEFLDGQRWIEDNLDVSVNSEVSVFEVNIRFIGGLLAAYYLSGEEIFKIKAVQLAEKLLP AFNTPGTIPWAMVNLKSGVGRNWGWASAGSSILAEFGTLHMEFIHLSYLTGDLTYKKVM HIRKLLQKMDRPNGLYPNYLNPRTGRWGQYHTSVGGLGDSFYEYLLKAWLMSDKTDHEAR KMYDDAIEAIEKHLIKSRGGLTFIGEWKNGHLEKKMGHLACFAGGMFALGADGSRADKA GHYLELGAEIARTCHESYDRTALKLPESFKFDGAVEAVAVRQAEKYIILRPEVIETYWY LWRFTDPRYRQWGWEAALAIEKYCRVNGGFGSGVKDVYSSTPTHDDVQQSFFLAETLKYL YLLFSGD DLLPLDHWVFNTEAHPLPVLHLANTTSLSGNPAVR
Gene Alias	MAN1B

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~73.00 kDa using MAN1A2 transfected 293T lysate. <u>Western Blotting (Transfected lysate)</u>
--------------------------	--



Western Blot analysis of MAN1A2 expression in transfected 293T cell line by MAN1A2 rabbit polyclonal antibody.

Lane 1: MAN1A2 transfected lysate(73.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619