

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-MAN2B1-248

Anti MAN2B1

Product type	Primary antibodies
Immunogen	MAN2B1 (NP_000519.2, 1 a.a. ~ 1011 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MGAYARASGVCARGCLDSAGPWTMSRALRPPLPLCFFLLLLAAAGARAGGYETCPTVQP NMLNVHLLPHTHDDVGWLKTVQYFYGIKNDIQHAGVQYILDSVISALLADPTRRFIYVE IAFFSRWWHQQTNATQEVVRDLVRQGRLEFANGGWVMNDEAATHYGAIVDQMTLGLRFLE DTFGNDGRPRVAWHIDPFGHSREQASLFAQMFGDGGFFGRLDYQDKWVRMQKLEMEQVWR ASTSLKPPTADLFTGVLPNGYNPPRNLCWDVLCVDQPLVEDPRSPEYNAKELVDYFLNVA TAQGRIYRTNHTVMTMGSDFYENANMWFKNLCLKLRLVNAQQAKGSSVHVLYSTPACYL WELNKANLTWSVKHDDFFPYADGPHQFWTGYFSSRPALKRYERLSYNFLQVCNQLEALVG LAANVGPGYSGDSAPLNEAMAVLQHHDVSGTSRQHVANDYARQLAAGWGCPCEVLLSNAL ARLRGFKDHFTFCQQLNISICPLSQTAARFQVIVYNPLGRKVNWMVRLPVSEGVFVVKDP NGRTVPSDVVIFPSSDSQAHPELLFSASLPALGFSTYSVAQVPRWKPQARAPQPIPRRS WSPALTIENEHIRATFDPTDGLLMEIMNMNQQLLPVRQTFFWYNASIGDNESDQASGAY IFRPNQQKPLPVSRWAQIHLVKTPLVQEVHQNFSAWCSQVVRLYPGQRHLELEWSVGPIIP VGDTWGKEVISRFDTPLETKGRFYTDSNGREILERRRDYRPTWKLNQTEPVAGNYYPVNT RIYITDGNMQLTVLTD RSQGGSSLRDGSLELMVHRRLKDDGRGVSEPLMENGSGAWVRG RHLVLLDTAQAAAAGHRLLAEQEVLAPQVVLAPGGGAAYNLGAPPRTQFSGLRRLDPPSV HLLTLASWGPPEMVLRLLEHQFAVGEDSGRNLAPVTLNLRDLFSTFTITRLQETTLVANQ LREAASRLKWTNTGPTPHQTPYQLDPANITLEPMEIRTF LASVQWKEVDG


Gene Alias LAMAN, MANB

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~113.70 kDa using MAN2B1 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of MAN2B1 expression in transfected 293T cell line by MAN2B1 rabbit polyclonal antibody.

Lane 1: MAN2B1 transfected lysate(113.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of MAN2B1 expression in HepG2 (Predicted Mw: 113.6730 kDa).

Western Blotting (Cell lysate)



MAN2B1 rabbit polyclonal antibody. Western Blot analysis of MAN2B1 expression in HepG2.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619