

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- MAP7-375

Anti MAP7


Product type	Primary antibodies
Immunogen	MAP7 (NP_003971.1, 1 a.a. ~ 749 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAELGAGGDGHRGGDGAVRSETAPDSYKVQDKKNASSRPASAI SGQNNNHSGNKPDP PPV LRVDDRQRLARERREEREKQLAAREIVWLEREERARQHYEKHLEERKKRLEEQRQKEERR RAAVEEKRRQRLEEDKERHEAVVRRTMERSQPKPKQKHNRWSWGGSLHGSPSIHSADPDRR SVSTMNLSKYVDPVISKRLSSSSATLLNSPDRARRLQLSPWESSVVRLLTPTHSFLARS KSTAALSGEAASCSPIMPYKAAHSRNSMDRPKLFVTPPEGSSRRRIHGTASYKKERER ENVLFLTSGTRRAVSPSNPKARQPARSRLWLPSKSLPHLPGTTPRPTSSLPPGSVKAAPAQ VRPPSPGNIRPVKREVKVEPEKKDPEKEPQKVANESLKGRAPLVKVEEATVEERTPAEP EVGPAAPAMAPAPASAPAPASAPAPAPVPTPAMVSAPSSTVNASASVKTSAGTTDPEEAT RLLAEKRLAREQREKEERERREQEELERQKREELAQRVAEERTTTRREEESRRLEAEQAR EKEEQLQRQAEERALREREEAERAQRQKEEEARVREEAERVQREREKHFQREEQERLERK KRLEEIMKRTRRTEATDKKTS DQRNGDIAGALTGGTEVSALPCTTNAPGNGKPVGSPHV VTSHQSKVTVESTPDLEKQPNENGVSQVQNFEEIINLPIGSKPSRLDVTNSESPEIPLN PILAFDDEGTLGPLPQVDGVQTQQTAEVI
	E-MAP-115, EMAP115

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~84.10 kDa using MAP7 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of MAP7 expression in transfected 293T cell line by MAP7 rabbit polyclonal antibody.

Lane 1: MAP7 transfected lysate(84.10 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619