

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP-MKAPK2-066

Anti MAPKAPK2

Product type	Primary antibodies
Immunogen	MAPKAPK2 (NP_116584.2, 1 a.a. ~ 400 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MLSNSQGQSPVPFPAPAPPPQPPTPALPHPPAQPPPPPQQFPQFHVKSGLQIKKNAII DDYKVTSQVLGLGINGKVLQIFNKRTQEKFALKMLQDCPKARREVELHWRASQCPHIVRI VDVYENLYAGRKCLLIVMECLDGGELFSRIQDRGDQAFTEREASEIMKSIGEAIQYLHSI NIAHRDVKPENLLYTSKRPNAILKLTDFGFAKETTSNLSLTTPCYTPYYVAPEVLGPEKY DKSCDMWSLGVIMYILLCGYPPFYSNHGLAISPGMKTRIRMGQYEFNPPEWSEVSEEVKM LIRNLLKTEPTQRMTEFMNHPWIMQSTKVPQTPLHTSRVLKEDKERWEDVKEEMTSAL ATMRVDYEQIKIKKIEDASNPLLLKRRKKARALEAAALAH

Gene Alias	MK2
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~45.60 kDa using MAPKAPK2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of MAPKAPK2 expression in transfected 293T cell line by MAPKAPK2 rabbit polyclonal antibody.

Lane 1: MAPKAPK2 transfected lysate(45.60 kDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of MAPKAPK2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of MAPKAPK2 transfected lysate using MAPKAPK2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with MAPKAPK2 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619