

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-MMP20-122

Anti MMP20

Product type	Primary antibodies
Immunogen	MMP20 (NP_004762.2, 1 a.a. ~ 483 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MKVLPAASGLAVFLIMALKFSTAAPSLVAASPRTWNNYRLAQAYLDKYYTNKEGHQIGEM VARGNSMIRKIKELQAFFGLQVTGKLDQTTMNVIKKPRCGVPDVANYRLFPGEKWKKN TLTYRISKYTPSMSSVEVDKAVEMALQAWSSAVPLSFVRINSGEADIMISFENGHDHGSY PFDGPRGTLAHAFAPGEGGLGGDTHFDNAEKWTMG TNGFNLFVAAHEFGHALGLAHSTDP SALMYPTYKYKNPYGFHLPKDDVKGIQALYGPRKVFLGKPTLPHAPHHKPSIPDLCDSSS SFDAVTMLGKELLFLKDRIFWRRQVHLRTGIRPSTITSSFPQLMSNVDAAYEVAERGTA FFKGPYHWITRGFQMGGPPRTIYDFGFPRHVQQIDAAYVLRPQKTLFFVGDEYYSYDER KRKMEKDYPKNTEEEFSGVNGQIDAAVELNGYIYFFSGPKTYKYDTEKEDVSVVKSSSW IGC
Gene Alias	MMP-20
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~53.13 kDa using MMP20 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of MMP20 expression in transfected 293T cell line by MMP20 rabbit polyclonal antibody.

Lane 1: MMP20 transfected lysate(53.13 KDa).

Lane 2: Non-transfected lysate.

I

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619