



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-MRPL13-563

Anti MRPL13

Product type	Primary antibodies
Immunogen	MRPL13 (NP_054797.2, 1 a.a. ~ 178 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSSFSRAPQQWATFARIWYLLDGKMQPPGKLAAMASIRLQGLHKPVYHALSDCGDHVVIM NTRHIAFSGNKWEQKVYSSHTGYPGGFRQVTAAQLHLRDPVAIVKLAIYGMLPKNLHRRT MMERLHLFPDEYIPEDILKNLVEELPQPRKIPKRLDEYTQEEIDAFRLWTPPEDYRL


Gene Alias L13, L13A, L13mt, RPL13, RPML13

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~20.70 kDa using MRPL13 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of MRPL13 expression in transfected 293T cell line by MRPL13 rabbit polyclonal antibody.

Lane 1: MRPL13 transfected lysate(20.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of MRPL13 expression in human kidney (Predicted Mw: 20.6920 kDa).

Western Blotting (Tissue lysate)



MRPL13 rabbit polyclonal antibody. Western Blot analysis of MRPL13 expression in human kidney.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619