

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-MTMR6-064

Anti MTMR6

Product type	Primary antibodies
Immunogen	MTMR6 (AAH40012.1, 1 a.a. ~ 621 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

**Immunogen
Sequence**

MEHIRTTKVEQVKLLDRFSTSNKSLTGTLTLTATHLLFIDSHQKETWILHHHIASVEKLA
LTTSGCPLVIQCKNFRTVHFIVPRERDCHDIYNSLLQLSKQAKYEDLYAFSYNPKQNDSE
RLQGWQLIDLAEYKRMGVPNSHWQLSDANRDYKICETYPRELYVPRIASKPIIVGSSKF
RSKGRFPVLSYYHQDKEAAICRCSQPLSGFSARCLEDEHLLQAISKANPVNRYMYVMDTR
PKLNAMANRAAGKGYENEDNYSNIRFQFVGIEINIHVMRSSLQKLLEVNGTKGLSVNDFYS
GLESSGWLRIKAVMDAAVFLAKAITVENASVLVHCSDGWDRTSQVCSLGSLLLDSEYRT
IKGFMVLIEKDWISFGHKFSERCGQLDGDPEKVPVFTQFLECWHLTEQFPQAFEFSEA
FLLQIHEHIHSCQFGNFLGNCQKEREELKLKEKTYSLWPFLLEDQKKYLNPLYSSSESHRF
TVLEPNTVSFNFKFWRNMYHQFDRTLHPRQSVFNIIMNMNEQNKQLEKDIKDLESKIKQR
KNKQTDGILTKELLHSVHPESPNLKTSLCFKEQTLLPVNDALRTIEGSSPADNRYSEYAE
EFSKSEPAVVSLEYGVARMTC


Gene Alias

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~72.00 kDa using MTMR6 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of MTMR6 expression in transfected 293T cell line by MTMR6 rabbit polyclonal antibody.

Lane 1: MTMR6 transfected lysate(72.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619