



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-MTP18-073

Anti MTP18

Product type Primary antibodies

Immunogen MTP18 (NP_057582.2, 1 a.a. ~ 166 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MSEPQPRGAERDLYRDTWVRYLGYANEVGEAFRSLVPAAVVWLSYGVASSYVLADAIDKG
KKAGEVPSPEAGRSARVTVAVVDTFVWQALASVAIPGFTINRVCAASLYVLGTATRWPLA
VRKWTTTALGLLTPIIIHPIDRSVDFLLDSSLRKLYPTVGKPSSS

Gene Alias HSPC242

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~18.00 kDa using MTP18 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of MTP18 expression in transfected 293T cell line by MTP18 rabbit polyclonal antibody.

Lane 1: MTP18 transfected lysate(18.00 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of MTP18 transfected lysate.

Immunoprecipitation



Immunoprecipitation of MTP18 transfected lysate using MTP18 rabbit polyclonal antibody and Protein A Magnetic Bead, and Immunoblotted with MTP18 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619