

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-NACA-182


Anti NACA

Product type	Primary antibodies
Immunogen	NACA (NP_005585.1, 1 a.a. ~ 215 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MPGEATETVPATEQELPQPQAETGSGTESDSDESVPLEEQDSTQATTQQAQLAAAAEID EETVSKAKQSRSEKKARKAMSKLGLRQVTGVTRVTIRKSKNILFVITKPDVYKSPASDTY IVFGEAKIEDLSQQAQLAAAEKFKVQGEAVSNIQENTQTPTVQEESEEEVDETGVKVD IELVMSQANVSRKAVRALKNNSNDIVNAIMELTM

Gene Alias HSD48, MGC117224, NACA1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~23.40 kDa using NACA transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of NACA expression in transfected 293T cell line by NACA rabbit polyclonal antibody.

Lane 1: NACA transfected lysate(23.40 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619